iHealth® View

Smart Blood Pressure Wrist Monitor

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USER GUIDE

Table of Contents

INTRODUCTION	1
PACKAGE CONTENTS	1
INTENDED USE	1
CONTRAINDICATION	2
PARTS AND DISPLAY INDICATORS	2
SET UP REQUIREMENTS	3
SET UP PROCEDURES	3
BODY POSTURE DURING MEASUREMENT	3
SYNC TIME AND RESULTS VIA <i>BLUETOOTH</i>	5
SPECIFICATIONS	6
GENERAL SAFETY AND PRECAUTIONS	7
BATTERY HANDLING AND USAGE	9
TROUBLESHOOTING	10
CARE AND MAINTENANCE	12
WARRANTY INFORMATION	13
EXPLANATION OF SYMBOLS	13
CONTACT AND CUSTOMER SERVICE	14
IMPORTANT INFORMATION REQUIRED BY THE FCC	14
OTHER STANDARDS AND COMPLIANCES	15

INTRODUCTION

Thank you for selecting the iHealth View Smart Blood Pressure Wrist Monitor. The iHealth View Smart Blood Pressure Wrist Monitor is a fully automatic wrist cuff blood pressure monitor that uses the oscillometric principle to measure your blood pressure and pulse rate. The monitor works with your mobile devices to track and share vital blood pressure data.

PACKAGE CONTENTS

- 1 iHealth View Smart Blood Pressure Wrist Monitor
- 1 User Guide
- 1 Quick Start Guide
- 1 Charging Cable
- 1 Travel Case

INTENDED USE



The iHealth View Smart Blood Pressure Wrist Monitor (Electronic Sphygmomanometer) is intended for use in a professional setting or at home and is a non-invasive blood pressure measurement system. It is designed to measure the systolic and diastolic blood pressures and pulse rate of an adult individual by using a technique in which an inflatable cuff is wrapped around the wrist. The measurement range of the standard cuff circumference is 5.3" to 8.7"(13.5cm-22cm).

Note: Consult your physician for proper interpretation of blood pressure results.

CONTRAINDICATION

⚠ It is not recommended for people with serious arrhythmia to use this Smart Blood Pressure Wrist Monitor.

PARTS AND DISPLAY INDICATORS

SET UP REQUIREMENTS

Compatible with iPhone 4s+, iPod Touch 5, iPad 3+, iPad Mini+, iPad Air+ and select Android phones. Requires iOS version 7.0+ and Android version 4.0+. For a complete list of compatible devices, visit our support page at www.ihealthlabs.eu/support.

SET UP PROCEDURES

Download the Free iHealth MyVitals App

Prior to first use, download and install “iHealth MyVitals” from the App Store or Google Play Store. Follow the on-screen instructions to register and set up your personal account.

Access the iHealth Cloud Account

Your iHealth account also gives you access to the free and secure iHealth cloud service. Go to www.ihealthlabs.eu and sign in with the same account.


Charge Battery before First Use

Connect the monitor to a USB port using the charging cable provided until the charging indicator  light green and  indicator steady.

BODY POSTURE DURING MEASUREMENT

Blood pressure can be affected by the position of the cuff and your physiologic condition. It is very important that the cuff is positioned at your heart level during blood pressure measurements.

1. Sit comfortably with your feet flat on the floor without crossing your legs. Stay still during measurement. Do not move your wrist, body, or the monitor.
2. Place your left hand palm-side up in front of you and leave 1-2cm between the monitor and the bottom of your left palm. If the monitor is correctly placed, iHealth logo will be facing upright. Please do not use on right wrist.

3. The center of the cuff should be at your heart level.

4. It is advised to place the monitor's travel case under your arm for support and to keep your arm at optimal height for measurement.

TAKING YOUR BLOOD PRESSURE READING


1. Apply the cuff or press the START/STOP button, the monitor will activate and begin to detect your wrist position.
2. Adjust the height of your wrist, the monitor will detect your wrist position and the measurement will start **ONLY** when the correct position is detected. When the position is correct, press "START/STOP" button again to start measurement.

Wrist position too low

Wrist position too high

Wrist position good

3. Then the cuff will be slowly inflated. The blood pressure and pulse will be measured during inflation. Inflation will stop as soon as the blood pressure and pulse rate have been calculated and displayed on the screen. The result will be automatically stored in the memory, and all results will be uploaded to the App automatically upon the next successful *Bluetooth* connection.

SYNC TIME AND RESULTS VIA *BLUETOOTH*

Connect to iOS Device via *Bluetooth*

1. Launch the iHealth MyVitals App from your iOS device.
2. Enable *Bluetooth* on your iOS device: Slide gently upward the iPhone screen from  to open the Control Center, and click on the .
3. Once the monitor is off, Press "START/STOP" button for 2 seconds until *Bluetooth* indicator begins flashing, when a successful connection has been established, the *Bluetooth* indicator light will stop flashing.
4. When sync is processing, the sync indicator flash  and  , when sync is finished, the sync indicator lights up .

Connect to Android Device via *Bluetooth*

- Press the START/STOP button for 2 seconds, the *Bluetooth* indicator will begin flashing.
- Enable the *Bluetooth* on your Android Device.
- When using the monitor for the first time, you should pair the monitor to the Android device. Go into your device's setting to pair it with your Android device. Check the *Bluetooth* menu for the model name of your monitor (BP7SXXX) to appear, and then select it to pair. This may take up to 30 seconds and your monitor's *Bluetooth* indicator will light up when a connection has been established.
- Launch the "iHealth MyVitals" app to sync with your monitor.
- When sync is processing, the sync indicator flash  and  , when sync is finished, the sync indicator lights up .

- Please repeat these steps when you switch to another Android device with the monitor.

Remember to:

1. Make sure your wrist size is within cuff circumference; see the range in “SPECIFICATIONS”
2. Measure on the same wrist each time.
3. Stay still and calm for one to one and half minutes before taking a blood pressure measurement. Prolonged over-inflation of the bladder may cause bruises of your wrist.
4. Keep the cuff clean. Cleaning the cuff after every 200 times of usage is recommended. If the cuff becomes dirty, clean it with a moistened cloth. Do not rinse the monitor or cuff with running water.
5. Press the “START/STOP” button at any time to interrupt a measurement.

Note: *Physical activity, eating, drinking, smoking, excitement, stress, and many other factors influence blood pressure results.*

SPECIFICATIONS

1. Product name: iHealth View Smart Blood Pressure Wrist Monitor
2. Model: BP7S
3. Classification: Internally powered; Type BF applied part; IP22, No AP or APG; Continuous operation
4. Machine size: approx. 2.8"×2.9"×0.7" (72mm×74mm×17.6mm)
5. Cuff circumference: 5.3"- 8.7" (13.5cm-22cm)
6. Weight: approx. 4.2oz(120g)(including cuff)
7. Power: DC:5.0V  1.0A, Battery: 1*3.7V  Li-ion 400mAh
8. Measurement range: Cuff pressure: 0-300 mmHg Systolic: 60-260 mmHg
Diastolic: 40-199 mmHg Pulse rate: 40-180 beats/minute
9. Wireless communication: *Bluetooth* V4.1 Frequency Band: 2.400-2.4835 GHz
10. Accuracy: Pressure: ±3 mmHg Pulse rate: ±5%

11. Environmental temperature for operation: 10°C-40°C(50°F -104°F)
12. Environmental humidity for operation: ≤85%RH
13. Environmental temperature for storage and transport:-20°C-55°C(-4°F-131°F)
14. Environmental humidity for storage and transport: ≤85%RH
15. Environmental pressure: 80kPa-105kPa
16. Battery life: more than 80 measurements on a full charge
17. The blood pressure measurement system includes accessories: pump, valve, cuff, and sensor.

Note: *These specifications are subject to change without notice.*

GENERAL SAFETY AND PRECAUTIONS

1. Read all of the information in the User Guide and other provided instructions before operating the unit.
2. Consult your physician for any of the following situations:
 - a) The application of the cuff over a wound or inflamed area.
 - b) The application of the cuff on any limb with intravascular access or therapy, or an arteriovenous (A-V) shunt.
 - c) The application of the cuff on the arm on the side of a mastectomy.
 - d) Simultaneous use with other medical monitoring equipment on the same limb.
 - e) The blood circulation of the user needs to be checked.
3. Do not use this product in a moving vehicle as this may result in inaccurate measurements.
4. Blood pressure measurements determined by this product are equivalent to those obtained by professional healthcare practitioners using the cuff/stethoscope auscultation method within the limits prescribed by the American National Standard, Electronic or Automated Sphygmomanometer. This device is also clinically validated according to the 2010 Protocol of the European Society of Hypertension (ESH 2010).
5. If an Irregular Heartbeat (IHB) is detected during the measurement procedure, the IHB symbol will be displayed in the “iHealth MyVitals” APP. Under this condition, the Smart Blood Pressure Wrist Monitor can keep functioning, but the results may be inaccurate. Please consult your physician for accurate assessment.

The IHB symbol will be displayed under 2 conditions:

(1) The coefficient of variation (CV) of pulse period >25%.

(2) The difference of adjacent pulse period is $\geq 0.14s$ and more than 53 percent of the total number of pulses readings falls within this definition.

6. Please do not use any cuff other than that supplied by the manufacturer as this may result in inaccurate measurements.

7. For information regarding potential electromagnetic or other interference between the blood pressure monitor and other devices together with advice regarding avoidance of such interference, please see ELECTROMAGNETIC COMPATIBILITY INFORMATION. It is suggested that the blood pressure monitor should be operated at least 10 meters away from electric or wireless devices (e.g. routers, microwaves oven, etc.)


8. If the blood pressure measurement (systolic or diastolic) is outside the rated range specified in part SPECIFICATIONS, the monitor will immediately display a technical alarm on the screen. In this case, repeat the measurement ensuring that the proper measurement procedures are followed and/or consult with your medical professional. The technical alarm is preset in the factory and cannot be adjusted or inactivated. This technical alarm is assigned as low priority according to IEC 60601-1-8. The technical alarm does not need to be reset.

9. This device requires a medical AC adapter with an output of DC 5.0V that complies with IEC 60601-1/UL 60601-1 and IEC 60601-1-2/EN 60601-1-2 such as ASP5-05010002JU (input: 100-240V, 50/60Hz, 200mA; output: DC 5V, 1.0A). Please note that the monitor jack size is USB mini B. The USB jack should be used for charging only.

⚠ This Monitor is designed for adults and should never be used on infants, young children, pregnant or pre-eclamptic patients. Consult your physician before use on children.

⚠ This product might not meet its performance specifications if stored or used outside the specified temperature and humidity ranges.

⚠ Please do not share the cuff with any infectious person to avoid cross-infection.

-  The patient is an intended operator.
-  Keep the device out of the reach of children/pets to avoid inhalation or swallowing of small parts.
-  If you are allergic to plastic/rubber, please don't use this device.

BATTERY HANDLING AND USAGE

- The battery charge will be displayed on the LED screen when the monitor is active. And when the monitor is connected to the “iHealth MyVitals” APP, the battery charge will be displayed in the APP. If the power is less than 25%, please charge the battery. The monitor will not work until the battery has enough power.
- When the monitor needs charging, please connect the monitor to a power source.
- You should charge the battery when the battery is less than 25% charged. Overcharging the battery may reduce its lifetime.
- When in charging mode, at first the monitor will activate and all display characters are shown for self-test. Please contact the service center if symbol is missing.

After the charging status will be displayed on the LED screen. See the table below for details.

Monitor Status	Status Indicator	
Charging	 symbol on	 symbol rolling
Fully charged	 symbol on	 symbol steady
Low battery	 symbol off	 symbol flashing (for a few seconds)

-  Do not change the battery. If the battery can no longer be charged, please contact Customer Service.
-  Overcharging the battery may reduce its lifetime.

- ⚠ Lithium battery replacement by inadequately trained personnel could result in a hazard such as a fire or explosion.
 - ⚠ Do not plug or unplug the power cord into the electrical outlet with wet hands. If the AC adapter is abnormal, please change the adapter.
 - ⚠ Do not pull out the adapter when you are using the monitor.
 - ⚠ Do not use any other type of AC adapter as it may harm the monitor.
 - ⚠ Do not use the monitor while charging.
 - ⚠ The monitor, cable, battery and cuff must be disposed of according to local regulations at the end of their usage.
- Note: Battery life and charge cycles vary by use and settings.*

TROUBLESHOOTING

PROBLEM	POSSIBLE CAUSE	SOLUTION
Low Battery	Battery do not have enough power	Charge the battery
LED display reads "Er0"	Pressure system is unstable before measurement	Retest, make sure not to move your arm or the monitor
LED display reads "Er1"	Fail to detect systolic pressure	
LED display reads "Er2"	Fail to detect diastolic pressure	
LED display reads "Er3"	Pneumatic system blocked or cuff is too tight during inflation	Apply the cuff correctly and try again
LED display reads "Er4"	Pneumatic system leakage or cuff is too loose during inflation	

LED display reads "Er5"	Cuff pressure above 300mmHg		
LED display reads "Er6"	More than 160 seconds with cuff pressure above 15 mmHg		
LED display reads "Er7"	memory accessing error		
LED display reads "Er8"	Device parameter checking error		
LED display reads "ErA"	Pressure sensor parameter error		
LED display reads "Er" 	<i>Bluetooth</i> communicate error		Connect the mobile device correctly and try again, if the monitor is still abnormal, please contact the local distributor or the factory.
The cuff position was not correct or it was not properly tightened	Review the cuff application instructions and retest		
LED display reads an abnormal result	Body posture was not correct during testing	Retest when calm; avoid speaking or movement during the test	
LED display reads an abnormal result	Speaking, moving arm or body, being angry, excited or nervous during test		
	<i>Bluetooth</i> connection unsuccessful, monitor is abnormal, or strong electromagnetic interference is present	Reset iOS/Android device. Reset monitor by pressing the START/STOP button and holding for about 10 seconds. Make sure the monitor and iOS/Android device are away from other electrical equipment. Please see GENERAL SAFETY AND PRECAUTIONS	
No response when you press button	Incorrect operation or strong electromagnetic interference	Press the START/STOP button and hold for about 10 seconds to reset the device.	

CARE AND MAINTENANCE

1. If this monitor is stored near freezing temperatures, allow it to return to room temperature before use.
2. If the monitor is not used for a long time, please be sure to fully charge it every month.
3. No monitor component needs to be maintained by the user. The circuit diagrams, component part lists, descriptions, calibration instructions, or other information which will assist the user's appropriately qualified technical personnel to repair those parts of the equipment which are designated for repair can be supplied by the iHealth technical department.
4. Clean the monitor with a dry, soft cloth or a moistened and well wrung soft cloth using water, diluted disinfectant alcohol, or diluted detergent.
5. The monitor can maintain the safety and performance characteristics for a minimum of 10,000 measurements or three years of usage, and the cuff integrity is maintained after 1,000 open close cycles.
6. The battery can maintain the performance characteristics for a minimum of 300 charge cycles.
7. It is recommended that if the cuff is used in a hospital or a clinic, it be disinfected twice a week. Wipe the inner side (the side that contacts skin) of the cuff with a soft cloth lightly moistened with Ethyl alcohol (75-90%). Then air dry the cuff.

- ⚠ Do not drop this monitor or subject it to strong impact.
- ⚠ Avoid high temperature and direct sunlight. Do not immerse the monitor in water as this will result in damage to the monitor.

⚠ Do not attempt to disassemble this monitor.

⚠ Not servicing/maintenance while the monitor is in use.

⚠ Battery replacement should only be performed by a qualified iHealth technician. To do otherwise will void your warranty and possibly damage your unit.

⚠ Cuff replacement should only be performed by a qualified iHealth technician. To do otherwise will possibly damage your unit.

8. It is recommended that product performance be checked every 2 years or after each repair. Please contact the iHealth Customer Service Center to do so.

WARRANTY INFORMATION

The iHealth View Smart Blood Pressure Wrist Monitor is warranted to be free from defects in materials and workmanship within one year from the date of purchase when used in accordance with the provided instructions. The warranty extends only to the end user. We will, at our option, repair or replace without charge the iHealth View Wireless Blood Pressure Wrist Monitor covered by the warranty. Repair or replacement is our only responsibility and your only remedy under the warranty.

EXPLANATION OF SYMBOLS

Symbol for "TYPE BF APPLIED PARTS" (cuff only)

Symbol for "MANUFACTURER"

Symbol for "THE OPERATION GUIDE MUST BE READ"

The sign background color: blue. The sign graphical symbol: white.

Symbol for "WARNING"

Symbol for "ENVIRONMENT PROTECTION – Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice".

Symbol for "SERIAL NUMBER"

Symbol for "EUROPEAN REPRESENTATIVE"

IP22

IP code of the device: this device's grade of ingress of solid foreign objects-->12.5mm diameter (and the against access to hazardous parts with finger); the grade of waterproof is dripping (15° tilted).

CE 0197

Symbol for "COMPLIES WITH MDD93/42/EEC REQUIREMENTS"
iHealth is a trademark of iHealth Labs, Inc.

iPad, iPhone, and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries.

CONTACT AND CUSTOMER SERVICE

Manufactured for iHealth Labs, Inc.
120 San Lucar Ct., Sunnyvale, CA 94086, USA
Tel: +1-855-816-7705 www.ihealthlabs.com

iHealthLabs Europe SAS www.ihealthlabs.eu
36 Rue de Ponthieu, 75008, Paris, France
E-mail: support@ihealthlabs.eu

ANDON HEALTH CO., LTD.
No. 3 Jinping Street, YaAn Road, Nankai District, Tianjin 300190, China.
Tel: 86-22-60526161

IMPORTANT INFORMATION REQUIRED BY THE FCC

This device complies with Part 15 of the FCC Rules. Its operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
 - (2) this device must accept any interference received, including interference that may cause undesired operation.
- Changes or modifications not expressly approved by iHealth Labs, Inc. would void the user's authority to operate the product.

Note: This product has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This product generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this product does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the

following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This product complies with Industry Canada. IC: RSS-210

IC NOTICE

This device complies with Industry Canada licence-exempt RSS standard(s). Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference that may cause undesired operation of the device.

Hereby, [ANDON HEALTH CO., LTD] declares that the radio equipment type [BP7S] is in compliance with Directive 2014/53/EU.

The full text of the EU declaration of conformity is available at the following internet address:www.ihealthlabs.eu

OTHER STANDARDS AND COMPLIANCES

The Smart Blood Pressure Wrist Monitor corresponds to the following standards:

IEC 60601-1Edition 3.1 2012-08/EN 60601-1:2006/A1:2013 (Medical electrical equipment – Part 1: General requirements for safety);

IEC60601-1-2:2014/EN 60601-1-2:2015 (Medical electrical equipment – Part 1: General requirements for safety; Collateral Standard-Electromagnetic compatibility - Requirements and tests);

EN 1060-1: 1995 + A1: 2002 + A2: 2009 (Non-invasive sphygmomanometers - Part 1: General requirements);

EN 1060-3: 1997 + A1: 2005 + A2: 2009 (Non-invasive sphygmomanometers - Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems);

IEC 80601-2-30 Edition 1.1 2013-07 (Medical electrical equipment –Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers).

ISO 81060-2 Second Edition 2013-05-01, Non-Invasive Sphygmomanometers - Part 2: Clinical Validation Of Automated Measurement Type.

ELECTROMAGNETIC COMPATIBILITY INFORMATION

This product is applicable to the equipment and system requirements for the purpose of receiving radio frequency energy for the purpose of the work, *Bluetooth* receive bandwidth 2M. This product can also be used to include RF transmitter equipment and system requirements and emission frequency of 2.4GHz ISM band, *Bluetooth* modulation types: GFSK , effective radiated power: $\leq 2.45\text{dBm}$

Table 1 - Emission

Phenomenon	Compliance	Electromagnetic environment
RF emissions	CISPR 11 Group 1, Class B	Home healthcare environment
Harmonic distortion	IEC 61000-3-2 Class A	Home healthcare environment
Voltage fluctuations and flicker	IEC 61000-3-3 Compliance	Home healthcare environment

Table 2 - Enclosure Port

Phenomenon	Basic EMC standard	Immunity test levels
		Home Healthcare Environment
Electrostatic Discharge	IEC 61000-4-2	±8 kV contact ±2kV, ±4kV, ±8kV, ±15kV air
Radiated RF EM field	IEC 61000-4-3	10V/m 80MHz-2.7GHz 80% AM at 1kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	Refer to table 3
Rated power frequency magnetic fields	IEC 61000-4-8	30A/m 50Hz or 60Hz

Table 3 – Proximity fields from RF wireless communications equipment

Test frequency (MHz)	Band (MHz)	Immunity test levels
		Professional healthcare facility environment
385	380-390	Pulse modulation 18Hz, 27V/m
450	430-470	FM, ±5kHz deviation, 1kHz sine, 28V/m
710	704-787	Pulse modulation 217Hz, 9V/m
745		
780		
810	800-960	Pulse modulation 18Hz, 28V/m
870		
930		
1720	1700-1990	Pulse modulation 217Hz, 28V/m
1845		
1970		
2450	2400-2570	Pulse modulation 217Hz, 28V/m

5240	5100-5800	Pulse modulation 217Hz, 9V/m
5500		
5785		

Table 4 – Input a.c. power Port

Phenomenon	Basic EMC standard	Immunity test levels
		Home Healthcare Environment
Electrical fast transients/burst	IEC 61000-4-4	±2 kV 100kHz repetition frequency
Surges Line-to-line	IEC 61000-4-5	±0.5 kV, ±1 kV
Surges Line-to-ground	IEC 61000-4-5	±0.5 kV, ±1 kV, ±2 kV
Conducted disturbances induced by RF fields	IEC 61000-4-6	3V, 0.15MHz-80MHz 6V in ISM and amateur radio bands between 0.15MHz and 80MHz 80%AM at 1kHz
Voltage dips	IEC 61000-4-11	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°
		0% UT; 1 cycle and 70% UT; 25/30 cycles Single phase: at 0°
Voltage interruptions	IEC 61000-4-11	0% UT; 250/300 cycles

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