
**IMPORTANT:** Refer to the iHealth Wireless Smart Gluco-Monitoring System's Owner's Manual for more detailed instructions on how to test blood glucose levels.

**PRE-SET UP INSTRUCTIONS**

- iOS device is compatible and is version V5.0 or higher.

**Wireless Smart Gluco-Monitoring System PACKAGE CONTENTS**


**DOWNLOAD THE FREE iHealth APP**

Prior to first use, download and install the free App “iHealth Gluco-Smart” from the App Store. Follow the onscreen instructions to register and set up your personal account.

**Connect the Glucose Meter to your iOS Device via Bluetooth**

- Press and hold the “Memory” button for 3 seconds to turn the meter on. At the same time, turn on *Bluetooth* function on your iOS device.
- In the *Bluetooth* menu, select the model name “BG5xxxxxx” and “Not Paired” to pair and connect. When the “” symbol appears on the meter, the meter is in *Bluetooth* communication mode.
- Each subsequent time you use the meter “Not Connected” will be displayed next to “BG5xxxxxx” in the *Bluetooth* Menu. Tap “Not Connected” to reestablish the *Bluetooth* connection.

**MEASUREMENT INSTRUCTIONS**


**Step 1: Open the Bluetooth capabilities of the IOS device**

**Step 2: Scan QR Code on Top of the Test Strip Vial**


When you are ready to open a new test strip vial, you must scan the QR Code on it. It is unnecessary to scan the QR Code for each test

**Step 3: Insert the test strip in the meter's Strip Port to turn on the meter**


**Step 4: Prepare the Lancing Device**


Pull off the lancet device cap

Insert a new lancet firmly into the lancet device cup holder


Twist off the lancet protective disk

Replace the lancet device cap

Set a lancet level

Lock the handle

Puncture your fingertip

### Step 5: Obtain a Blood Sample

Press the lancing device against the site to be lanced. Press the release button to puncture your fingertip. Gently squeeze and/or massage your fingertip until a round drop of blood at least 0.7 micro liter forms on your fingertip.

**Note:** For Alternative Site Testing and Control Solution Testing, it is recommended that you read the iHealth Wireless Smart Gluco-Monitoring System's Owner's Manual.

### Step 6: Apply the Blood Sample to the Test Strip

Apply the blood sample to the absorbent hole of the test strip. Make sure the confirmation window of the test strip is completely filled with the blood sample. Follow the on-screen instructions in the App to complete the blood glucose testing.

### Step 7: Remove Test Strip from the Blood Sample

Your iOS device will beep when you should remove the test strip from the blood sample.


### Step 8: Read the Test Results

Wait for about 5 seconds, your blood glucose result will appear on the iOS device screen.

### Step 9: Discard the Used Test Strip and Lancet

Remove the used test strip from the glucose meter using a small amount of tissue paper. Discard the used test strip and lancet properly. (Tip: Stick the lancet in the protective disk and dispose properly.)

## Taking offline Measurements without an iOS Device (Meter only)

After first successful *Bluetooth* connection, you may perform a test even when the meter is not connected to your iOS device via *Bluetooth*, the meter continuous offline measurement 25 times and the meter will save up to 500 test results.

iHealth is a trademark of iHealth Lab Inc.  
©2010 iHealth Lab Inc. All rights reserved.

Manufactured for iHealth Lab Inc.

USA:

iHealth Lab Inc. [www.ihealthlabs.com](http://www.ihealthlabs.com)

Mountain View, CA 94043, USA

+1-855-816-7705 (8:30 AM-5:00 PM PST, Monday to Friday)

Europe:

iHealthLabs Europe SARL [www.ihealthlabs.eu](http://www.ihealthlabs.eu)

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iHealth ist eine Marke der iHealth Lab Inc.
©2010 iHealth Lab Inc. Alle Rechte vorbehalten.
Hergestellt für iHealth Lab Inc.