▸ Welcome

Congratulations on your new iHealth Wireless Scale Lina with *Bluetooth*® Low Energy technology. This Scale tracks your weight and BMI, along with the free companion “iHealth MyVitals” app, can help you stay motivated and reach your weight goals.

▸ Mobile Device Compatibility

Works with the following devices:

iPhone 4S+

iPad 3+

iPad Mini+

iPad Air+

iPod Touch (5th generation)

Select Android devices

Requires iOS version 9.0+ or Android version 4.4+

Compatible devices are subject to change. For the latest compatibility list, (for US) go to www.ihealthlabs.com/support or <http://www.ihealthlabs.eu> (for EU)

▸ Companion Mobile App

Scan the QR code to download the free “iHealth MyVitals” companion app or directly from the Apple App Store or the Google play. Follow the on-screen instructions to register and set up your iHealth ID.

iHealth MyVitals

▸ Access iHealth Cloud

Upon setting up your iHealth ID, you will also have access to a free, secure iHealth cloud account. Go to www.ihealthlabs.com (for US) or <http://www.ihealthlabs.eu> (for EU), then click on “Sign In” to access your cloud account from a PC or Mac using your iHealth ID.

▸ Install the Batteries

Open the battery compartment on the bottom of the Scale and install the four batteries. The battery power should last approximately 3 months when measuring two times a day.

▸ First Time Use


For your first weigh-in, follow these steps to ensure the Scale and the app are connected and that the date and time synchronize correctly. If you have an iOS device:

1. Enable the *Bluetooth* on your mobile device (Settings->*Bluetooth*->On).

2. Tap the Scale to turn it on.

3. Launch the app to initiate the connection. The *Bluetooth* icon on your device will light up and stop flashing when a successful connection is established. You will also see the *Bluetooth* icon on your Scale lights up. The Scale is now connected and ready for your weigh-in.

• **Your Next Weigh-in**

Subsequent daily weigh-ins can be taken without being connected to your mobile device, as your weight data is stored in the Scale and uploaded to the app upon your next connection. Real-time weight results are viewable in the app if the Scale is connected to your mobile device and the app is running when your weigh-in occurs.

• **Software updates**

From time to time, the Scale's firmware will be updated to enhance performance. When a software update is available, a notification message will appear in the app. Complete the following steps to install the update:

- ➊ Follow steps in the "First-time Use" section.
- ➋ You will see a notification message in the app asking for permission to proceed. Select "Yes".
- ➌ The update will begin to download. You will see a cursor on the scale display moving from left to right. The scale will display "0.0" when the update is complete.

Visit <http://www.ihealthlabs.com> (for US) or <http://www.ihealthlabs.eu> (for EU) to obtain additional product information.

For Customer Service, please call +1(855) 816-7705 (USA).

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